



Alltesta™ Mini-Autosampler

Technical Specifications

Technical Specifications

Parameters	Alltesta™ Mini-Autosampler (Standard configuration)	Alltesta™ Mini-Autosampler (Customizable configuration)
Dimensions (WHL)	160 x 145 x 195 mm / 6.3 x 5.71 x 7.68 in	depends on configuration
Weight	3 kg / 6.6 lb	depends on configuration
Power	24 V	24 V
Plate	48 vials	48 vials 96-well plate 384-well plate
Pressure	Up to 5000 psi (345 bar)	depends on configuration
Syringe Capacity	4000 µL low-pressure stainless	120 µL high-pressure 4000 µL low-pressure stainless 4000 µL low-pressure metal free
Valve options	6x2	6x2 7x6 10x9
Needle	Short (36mm)	Short (36mm) Long (46mm) Capillary metal free (36mm)
Volume Accuracy	0.1 µL	depends on configuration
Communication	RS232 Virtual COM port via USB	RS232 Virtual COM port via USB