



# **Alltesta™ Mini UV-Vis Detector**

## **Technical Specifications**

## Technical Specifications

<b>Parameters</b>	<b>Alltesta™ Mini UV-Vis Detector (Standard configuration)</b>
Dimensions (WHL)	63 x 46 x 116 mm / 2.48 x 1.81 x 4.57 in
Weight	0.25 kg / 0.55 lb
Optical Cell Volume	Available in <3 µL and 10 µL options
Cell Optical Path	7.0 mm (other types of cells are available)
Liquid Ports	10-32 UNF compression
Built-in pressure regulator	Adjustable, 10-75 psi internal pressure
Maximum cell pressure	250 psi (18 bar)
Wavelength Range	235 nm and up
Simultaneous Wavelengths	Up to 4
Data Rate	20 sps/channel
Noise Level	0.03 mAU
Analog signal output	- 1.0 V to + 1.0 V
Power	USB Type B